FLOCHECK F – Instruktion

Dragbänksvakt

Om utrustningen

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
Down-flow Bench safety guard

About the equipment

FLOCHECK F is used to set the bench in normal or standby mode. The panel should also alert personnel to an incorrect airflow. Connected controller maintains the airflow that ensures that environmentally harmful substances are evacuated with exhaust air without leakage to the room. Setpoints and alarm limits are set according to current regulations. Changes to settings may only be made by a person skilled in the art. In case of suspicion of malfunction, contact a safety representative or a supervisor. Deviations from these specified values of settings may occur.

Function

FLOCHECK F shows the current status and mode of operation. The unit alerts by text messages, summer tone and by the activation of LEDs. The rate of airflow is set manually by pressing <ON> for the normal workflow or <OFF> for the stand-by rate of airflow.

Operating mode

Normal flow (Working mode)

Green LED and the text "OK" in the display's status bar indicates that the airflow is within the specified limits. The bench is activated to the required workflow. Operating mode "NORM" appears in the status bar.

Flashing red LED *with* summer tone warns that the airflow is lower than the intended workflow. The summer can temporarily be switched off by pressing <Mute>. As a reminder, the buzzer returns after 20 minutes.

Flashing red LED *without* summer tone indicates a higher workflow than intended.

When not in use: Switch to *standby mode* by pressing <OFF>. The down-flow bench maintains it's standby rate of airflow while reducing energy usage and wear.

Minimum flow (Standby mode)

Solid red LED alerts the operator that the down-flow bench is not in its working mode. The airflow is reduced to and maintained at the lower standby level. The operating mode "MIN" is displayed in the status bar. Alarms when the stand-by airflow is higher or lower than intended is indicated only by flashing red LED.

Settings

Settings can be changed by an authorised technician using the code protected menu system in FLOCHECK F.